

BioPAK Version 8

DIGITAL DENTAL

UPGRADE TO THE LATEST VERSION OF BIOPAK



Available to any current user of BioPAK Software - operating JVA, JT and EMG and network versions.

Silver support package = one year of free upgrades, rental units available for recalibrations and repairs

If you are upgrading to a 64 bit computing you need this upgrade

Limited one off price

- Upgrade from any previous version for \$950 and receive a 12 month Silver Support Package.
- Normally the cost to upgrade is \$950 per version (eg from Version 6 to Version 8 = \$1,900)
- Offer expires 28 February 2014

Major improvements include

Windows 7 and 8 compatible
Only version compatible with 64 bit computers
JVA "Doppler simulation" included
Features online continuous program updates over the year

OTHER CHANGES

- Added support for *Dentrix* software link.
- Added a new mastication ACP Sweep view that allows time synchronization with the ACC View.
- Added the JVA magnification to the JVA Sweep graph when the Grid is displayed.
- Improved CO point detection for EGN traces.
- The X-Y cursor now moves when the EGN Sweep view cursor is moved.
- Changed clench onset timing measurements in the EMG Summary view to display zero for the first muscle instead of a dash. Dashes are still used for muscles that don't exceed the threshold.
- Added a new view to EMG mastication analysis that displays separate RMS graphs for each muscle, showing the selected cycles drawn on top of each other.
- Added support for the Measurement Computing USB-1608G 16 channel data acquisition board, allowing recording of jaw tracking with up to 8 muscles simultaneously.
- Migrated to Firebird database engine version 2.5.
- BioPAK displays the data acquisition hardware detected after entering record mode.
- Improved the Acceleration and Jerk calculations in the ASCII conversion feature.
- Note: T Scan and BioTENS are not part of the BioPAK software system

VISIT OUR
WEBSITE



DOWNLOAD
MORE INFO



CALL US
07 3103 3373



EMAIL US IF
INTERESTED

